



CYP2C19

# Genedrive® CYP2C19 Collection Kit

The **Genedrive® Collection Kit** enables patient samples to be collected from the patient and **tested where required**.

## Genedrive® Collection Kit Procedure



Firmly swab  
the patient



Vigorously rotate the  
swab to release the DNA



Complete patient ID label  
attach to collection tube

## Send for testing



Tested on the  
Genedrive® System

**\*2/\*3  
Poor  
Metaboliser**

Results to inform  
treatment decisions

It is intended to be used  
only in combination with the  
**Genedrive® CYP2C19 ID Kit**  
and the Genedrive® System.




For more information  
scan the QR code

The **Genedrive® CYP2C19  
Collection** kit includes:



1 x Lysis buffer



1 x Patient ID label



1 x Buccal Swab

## Intellectual Property

\*References available on request. Genedrive® is a registered trademark of Genedrive Diagnostics Ltd.  
Copyright © 2025 genedrive plc. All rights reserved. GPM-187 V1 October 2025  
For further information on CYP2C19 genetic testing please refer to GPM-136

**genedrive**